



[www.nausicaa-medical.com](http://www.nausicaa-medical.com)

This medical device is intended for patient positioning. This is not a restraint system.

Before each use, it is important to check: the presence of the label, the stitches of the straps: that they do not fray or break, that the device does not show tears, snags, holes or signs of fraying (fabric and straps), that the buckles are not cracked or have no weak fastening.

## LARGE ABDOMINAL BELT

FOR ALL CHAIRS



TO COMPENSATE SLACKENING  
OF ABDOMINAL MUSCLES TO  
POSITION IN THE BED



Maximum useful area: + 230 cm  
Height of the belt: + 18 cm  
Buckle: 50 mm

Its large contact surface prevents any abdominal constriction.

Available in 4 sizes S, M, L, XL

Maximum useful area: 230 cm


Waist (in cm)



The use of the device requires the presence and vigilance of a third person. When used by an agitated or disoriented person, monitoring becomes imperative.



This medical device must be handled only by employees trained (by a healthcare professional) and qualified for its use (intended purpose, installation and environment of use, etc.). If the medical device is incorrectly installed, there is a risk of accident for the patient.



CCA2HP-PF

v 07.02/25

### Important Recommendations :

- In order to get a maximum efficiency of these sort of product, it is necessary :
- to choose the appropriate size for the patient
- to get the best possible adjustment to the patient
- These products must not be in direct contact with a wounded skin.



## LARGE ABDOMINAL BELT

### HOSPITAL RANGE

Combination of polyurethan  
foam and on both faces  
Jersey/PU coating  
100% polyester bias  
100% polyester strap  
Nylon buckle

- Place the belt around the patient's waist, smooth face outwards.
- Put the pointed end through the slit.
- Pass the strap around the backrest of the chair or around the mattress.
- Adjust the length of the straps through the Nylon buckles.
- Used in bed, the degree of restraint varies depending on the setting :
  - Slit on the stomach : kept on the back without any possibility for the patient to reach the buckles.
  - Slit in the back : possibility to move laterally.