



Before each use, it's important to check: that the label is attached; that the straps are not frayed or coming undone; that the device is free of tears, snags, holes or signs of fraying (fabric and fastening straps); and that the buckles are not cracked or showing signs of weak attachment.

RAISING SLING



Model Compositions:
Sterithane laminate
40 mm polyester webbing –
red 22 mm
Polyester bias
- black
Polyamide handle

This technical aid is designed for caregivers to help them lift a patient with a low level of physical activity into bed with less effort and without damaging their backs.

Product benefits :

- Helps combat RSI linked to carrying heavy loads and/or poor posture.
- Rigid handles for greater comfort and a better grip.
- Easy to use and clean



Before each use, it's important to check: that the label is attached; that strap stitching is not loose or frayed; that the device is free of tears, cracks, holes or signs of fraying (fabric and hooking straps); and that buckles are not cracked or showing signs of weak attachment.



Use of the device requires the presence and vigilance of a third party. In the event of use by an agitated or dis-orient person, supervision becomes imperative.

This medical device must be operated or trained by a health-care professional.



SRL-PF



v 01.05/24



Recommendations for use:

Single-person use:

The raising sling should be hooked to one side of the bed at a fixed point towards the top of the bed (slatted frame, gallows). The padded sleeve should be positioned at the patient's buttocks. The caregiver retrieves the other end of the sling from the other side of the bed and uses a splitting movement to raise the patient. If able, the patient can participate by looking at his stomach and pulling on the gallows handle and/or pushing on his feet.

Two-person use:

The sheath is placed on the patient's buttocks, with one person on each side of the bed holding the strap by its handles. The two people synchronize their movements to lift the patient to the top of the bed. As mentioned in the previous paragraph, the patient can participate to the extent of his or her ability.

Use with sliding sheet :

The raising sling can be used in combination with a sliding sheet, which is fitted beforehand. This combination enables effortless raising of heavier patients or those with less physical capacity to participate.