# NAUSICAA

## MEDICAL

www.nausicaa-medical.com

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EN

USER MANUAL  
PUMP: **R'GO SOINS** ——— p 24-46

© Nausicaa Médical S.A.  
Design Policy and Copyright:

As our policy is constantly evolving, we reserve the right to change designs without notice. It is forbidden to copy the content of this publication in whole or in part without the express permission of NAUSICAA MEDICAL.

# Summary

## User manuel / Pump : R'GO SOINS

1. Introduction	26
2. Content of the box	26
3. Safety information and instructions	27
3.1 Symbols	27
3.2 Important	27
3.3 Safety instructions	28
3.4 Precautions	29
3.5 Electromagnetic Compatibility (EMC)	29
4. Target groups	29
4.1 Targets patients	29
4.2 Users	29
5. Operating principle	30
6. Clinical Principles	30
6.1 Objective	30
6.2 Expected clinical benefit	30
6.3 Indications	30
6.4 Contraindications	30
6.5 Other measures required	31
6.6 Effects	31
6.7 Restrictions	31
6.8 Manual control	31
7. Presentation	32
7.1 The adjustment buttons	32
7.2 Description	32
8. Installation	33
8.1 Procedure	33
9. Getting started	33
9.1 Precautions and Restriction of Use	34
10. Settings	34
10.1 Comfort adjustment	34
10.2 Cycle time adjustment	35
10.3 Mode setting	35
10.3.1 Dynamic Mode	35
10.3.2 Static Mode	35
10.3.3 Care Mode	36
10.3.4 Sitting position	36
11. Uninstalling	37

# Summary

---

12. Alarms	37
12.1 LOW PRESSURE	37
12.2 SERVICE	37
12.3 SYSTEM FAILURE	38
13. Maintenance	38
13.1 Daily or weekly maintenance	38
13.1.1 Procedure	38
14. Pre-storage maintenance	39
14.1 Generality	39
14.2 Procedure	39
15. Storage	40
16. Climatic requirements	40
17. Compatible supports	41
18. Troubleshooting	42
19. Symbols (product label and/or packaging)	43
20. Technical data	44
21. Elimination	44
22. Lifespan	44
23. Warranty	45

# Introduction

## 1. Introduction

Dear Customer,

Nausicaa Medical thanks you for choosing this pump.

 To avoid any damage or hazards caused by misuse, we invite you to read this manual in its entirety before first use and before any subsequent use. This manual contains important information and remarks necessary for the proper use of this compressor. If you have any questions about safety warnings, precautions, or use, contact your local service provider or Nausicaa Medical's after-sales service. To avoid injury and damage due to misuse, do not use the compressor until all your questions have been answered.

 All users must have read this manual in its entirety before using this pump.

The service life of this pump depends mainly on proper installation and proper use.

It meets the following norms:

EN 60601-1 : 2006/A1 : 2013

EN 60601-1-2 : 2015

EN 60601-1-11 : 2015

EU Declaration of Conformity 2017/745

 The term “support” in this manual refers to the overlay or mattress used.  
The term “system” refers to all the elements that allow it to be used, i.e.: “this pump and a support”.


## Content of the box

### 2. Content of the box

Please check that the content listed below corresponds to the content you received.

1 digitally adjustable pump

1 user's manual 

# Safety information and instructions

## 3. Safety information and instructions

### 3.1 Symbols

To make this user manual easier to read and understand, here is a brief explanation of the main symbols used. Some of these symbols and their explanations refer to special attention.

Indicates that a manual is available.

Indicates that the manual should be read before use.

INFO :

Provides advice or information that the user and/or patient should be aware of.

**WARNING/DANGER:**

**Indicates an appropriate use, operation, or maintenance procedure to prevent tampering with and/or rendering inoperative this pump or any of its components or damage to other property.**

**Attracts the attention to a potential hazard requiring appropriate procedures or practices to prevent injury.**

### 3.2 Important

Read and follow the safety instructions carefully. They consist of a text or a combination of a symbol and a text. **The symbol used does not replace the text!**

For more information on the symbols used, see page 43.

The safety instructions apply to any person using and/or intervening in any way whatsoever on this pump or any of its components.

The person making this pump available to the user and/or a patient must inform the user of the potential risks and train them in the correct use so that no additional risks result from inappropriate behaviour.

**This pump should never be used unattended, especially when connected to the mains.**

**Be vigilant when children, people with disabilities or pets are present near the pump.**

**The presence of a third person is strongly recommended, especially if the pump triggers an alarm.**

**Use this pump only for its intended purpose. Any other use is prohibited.**

**DO NOT use with any other medium than those compatible and marketed by Nausicaa Medical.**

# Safety instructions

## 3.3 Safety instructions

**⚠ The electrical equipment (cable and plug) of this pump can be dangerous if improperly handled.**

The plug must always be accessible so that, if necessary, it is possible to completely disconnect the pump from the mains.

Never use the appliance if the plug and/or cable are damaged.

Do not use the pump if it has been dropped on the floor or in liquid, if it is damaged, if it has significant wear, if it makes an unusual noise, if it does not work properly, or if you are unsure of its safety. Always unplug the pump immediately after use.

Repairs and/or interventions on this pump must be carried out only by qualified personnel or by personnel who have received the training necessary. No modifications are allowed, and only original parts should be used.

To ensure the safety of the pump, the electrical components have been tested by an external and independent testing institute. However, improper use can still lead to the risk of a short circuit or electric shock.

If the pump has fallen into a liquid, never retrieve it directly without first carefully removing the plug from the wall outlet connected to the mains.

Do not use the pump near flammable liquids or gases. Protect it from direct sunlight. Keep it away from all sources of heat (heating, electric blanket,..) and open flames (candles,...).

Note that the pump draws in the surrounding air, so any smoke and/or excess dust can interfere with the degradation of its internal components.

To prevent the risk of overheating, never insert or use the pump in an enclosed compartment (box, etc.), never cover it with anything (sheet, blanket, etc.) and never position it on a soft and unstable surface to prevent the suction holes from always becoming clogged.

To prevent damage to the electrical cable, avoid all moving mechanisms (bed rails, frames, etc.).

If the installed support has an electrical cable management, it must be used.

Do not use any components other than those recommended by Nausicaa Medical.

The use of a non-compliant component that cannot guarantee the proper functioning of the pump systematically results in the cancellation of the legal warranty.

**i** It is recommended that you report to the service provider, the manufacturer, and the competent authority any serious incident\* that has occurred during the use of this pump or if you consider or have reason to believe that this device has been falsified.

\*»serious incident» means any incident that has directly or indirectly led to, or could lead to:

(a) The death of a patient, user, or other person.

(b) To the serious deterioration, temporary or permanent, of the state of health of a patient, user or other person.

# Safety instructions

---

## 3.4 Precautions

Protect the pump from sharp objects and surfaces, including animal claws or teeth.

Before use or reuse, the pump should be cleaned, disinfected, checked, clean and dry. Observe the recommendations for maintenance, cleaning, and disinfection (pages 38-40).

-  Store the pump as directed in this user manual (page 40).
-  Respect the minimum or maximum weight supported by the pump. Failure to comply with these weights can be very harmful to the patient. In addition, exceeding the maximum weight allowed can cause malfunction and internal damage to the pump, making it inoperable.

## 3.5 Electromagnetic Compatibility (EMC)

This pump complies with the standards of electromagnetic compatibility of electromedical devices.

Please note that possible electromagnetic or other influences between the pump and other electrical devices cannot be completely excluded.

If there is a risk of mutual interference, remove or move away (min. 3 m) from the source of interference.

Contact NAUSICAA MEDICAL for more detailed information about EMC.

# Target groups

---

## 4. Target groups

### 4.1 Target patients

The target patients are mainly adults with a height of > to 150 cm and a weight between 30 kg minimum and 200\* kg maximum and who, due to illness, injury, disability or age, have a bed rest period of more than 15 hours per day, have one or more pressure ulcers or are considered to be at risk due to the temporary or permanent deterioration of their state of health.

\*Depending on the support used

### 4.2 Users

Users are generally health professionals, assisted by caregivers such as: caregiver, carer or family member.

That's why this pump has been designed to make it easy and convenient to use and maintain, regardless of the user.

Regarding use, it is essential that all users must have received all the appropriate instructions for the proper use of this pump and that they can perform the supervisory function.

# Operating principle

## 5. Operating principle

The operating principle is the alternation of inflation of the cells of the support used according to a specific pattern, thus creating a cyclic pressure movement, and alternating on the different areas of the body.

## Clinical Principles

ENGLISH

## 6. Clinical Principles

The pump is one of the components of a **system** intended to be used for medical purposes to prevent the appearance of a pressure ulcer and/or to assist in its care.

### 6.1 Objective

Reduce pressure on the areas of the body that are most prone to the development of a pressure ulcer.

### 6.2 Expected clinical benefit

Maintenance of tissue oxygenation in the anatomical areas in contact with the support by reducing the pressure applied to the skin and subcutaneous tissues.

### 6.3 Indications

This **system** is mainly recommended in the context of the prevention and/or treatment of all types of bedsores (according to medical advice).

It is intended to reduce the incidence of pressure ulcers while optimizing patient comfort. It is an appropriate solution for patients who are bedridden for more than 15 hours a day and up during the day whose risk assessment is classified as «medium to high\*». As an aid to care, it is indicated for a constituted pressure ulcer judged to be stage 4\*.

All these indications are part of a clinical judgment as well as an individual protocol including repositioning, nutritional support, skin care, wound care, etc. All other aspects of care should be considered by a doctor as well as a senior caregiver.

\*Following medical advice, evaluation of the score via the Braden scale and depending on the support used.

**i** Regularly check the adequacy of the **system** with the patient's health status, including whether the patient's physical condition has changed (amputation, weight gain/loss, mental changes, etc.).

Always assess the patient's health status with a holistic view.

Always consider contraindications. If necessary, adjust the pump settings or choose a **more suitable** system.

### 6.4 Contraindications

Do not use the **system** if the patient has one or more of the following contraindications:

- Less than 30 kg or more than 200 kg\* (\*depending on the support used)
- Unstable spinal injury or other spinal disorder,
- Cervical or skeletal traction,
- Unstable spinal cord injury,
- Acute multiple traumas,
- Unstable posttraumatic bone fracture.

# Clinical Principles

If this is the case, use is subject to the advice of a qualified doctor.

- i** It is also recommended that this system not be used in a hyperbaric chamber or on a stretcher.

## 6.5 Other measures required

The **system** alone is not enough; other preventive measures are also essential, including:

- Changing position (at least every 2 to 3 hours);
- Maintain good skin hygiene and avoid or reduce maceration as much as possible;
- Change sanitary products regularly, especially in cases of severe incontinence;
- Observe or have observed the skin condition on a daily basis;
- Check that the food is sufficient and appropriate;
- Check that the patient is hydrating regularly and in sufficient quantities.

## 6.6 Effects

There are no scientifically substantiated adverse effects, although there may be a small risk that the alternating action of the pump will cause and/or promote spasticity. In these cases, the opinion of a specialist doctor should be sought.


## 6.7 Restrictions

- !** The pump is only suitable for people whose weight and height are within the maximum permissible limit (page 22 technical data). It is the different settings offered by the pump that determine the maximum weight supported.

- i** All the above (indications, contraindications, precautions, restrictions, effects, etc.) is only a guideline and can never replace a doctor's judgment or experience!

## 6.8 Manual control

In alternating mode, you can check if the pressure is adjusted correctly. To do this, locate the least swollen cell under the patient's buttocks and try to insert a flat hand with the palm facing up (see illustration below).

You should be able to feel a slight pressure while still having the ability to move back and forth. The acceptable range between the two is about 25 to 40 mm.

This manual verification procedure is issued by the AHCPR (Agency for Health Care Policy and Research).

- i** If you are completely unable to perform the back-and-forth movement, it means that the pressure is too low (bottoming out) > **increase the pressure.**
- i** If you can perform this movement without any resistance or even with extreme ease, it means that the pressure is too high > **reduce the pressure.**

# Presentation


## 7. Presentation

The pump can be used in both in-house and extramural acute and long-term care settings, in nursing homes, etc. as well as at home.

Its control panel is digitally adjustable. Its electronic management gives it, among other things, automatic pressure maintenance.

### 7.1 The adjustment buttons

1. On/off switch
2. Comfort pressure adjustment
3. Cycle time adjustment
4. Resetting alarms
5. Static or dynamic mode control
6. Activating Care Mode
7. Control-panel locking/unlocking

### 7.2 Description

- **On/Off switch (1) :**

Starts or stops the pump operation.

- **Comfort pressure adjustment (2) :**

Eight settings to optimize patient comfort, including 1 sitting position.

- **Cycle time adjustment (3) :**

Four values available: 10, 15, 20 and 25 minutes.

- **Alarms (4) :**

Puts the audible alarm or visual alarm to standby.  
Resets alarms.

- **Static or Dynamic Mode (5) :**

Enables/disables the selected mode.

- **Care Mode (6) :**

Enables/disables this mode.

**i** For safety, this mode automatically cancels after 30 minutes and returns to the previously selected setting.

- **Locking / Unlocking (7) :**

Manually locks/unlocks the control panel.

**i** The control panel locks automatically after 5 minutes of non-use. The chosen settings block and no changes can be made without manual intervention. To manually lock or unlock the control panel, press the button (7) for  $\pm 3$  seconds, and then release.

# Installation

## 8. Installation

Before installation, always check the pump for damage.

**DO NOT install or use a pump, even slightly damaged!**

### 8.1 Procedure

- Using the hanging hooks, hang the pump on the panel of the footboard.

- i** If this is not possible, hang it from the frame, making sure that it does not cause vibration, BUT NEVER place it on the ground.
- Place the power cord so that it does not pose any danger, so that it cannot be caught, crushed, jammed, or torn.
- i** If the support used is equipped with a cable duct, it is strongly recommended to insert the power cable into it.
- **Insert the plug of the power cable into the wall outlet.**
- i** Make sure that the electrical specifications of the sector are in line with those necessary for the operation of the pump.

## Getting starting

### 9. Getting started

**Step 1** ► Turn switch (1) to the «I» position. The green light comes on and a brief beep is heard. First, the pump enters the power-up control mode, and its software resets its settings. It is only after performing these actions that it begins to release air.

- i** When it is started, the pump uses static mode until it has blown 15 mmHg of air into the support used and then it automatically switches to alternating mode.

**Step 2** ► Allow the air to escape from the pump for a few seconds, then insert the 3-way connector of the air supply hoses from the support into the connector of the pump and let the support reach the comfort pressure «**STANDARD**> 5th LED» (page 36).

# Getting started

**⚠ Make sure that both connectors (support and pump) are properly connected and secured.**

**i** Be aware that the pressure application time depends on the support used  
Example: it takes about **25** minutes for an 88x13x200 cm support to be fully operational.

**i** **Using the Nausicaa Rapid Inflation Pump allows you to reach the required pressure in ± 2 minutes.**

**⚠ If the pressure of 15 mmHg is not reached within 70 minutes, a visual and audible alarm is triggered to warn of the problem (pages 39 & 43).**

**Step 3** ▶ Once the pressure is reached, the **system** is ready for use. You can then place the patient on the support and adjust the comfort pressure (see 10.1) if necessary.

## 9.1 Precautions & Restriction of Use

Take all necessary precautions

- When the pump is used at home without the intervention of nursing or similar staff (nursing assistant, etc.).
- When used in the first few days of post-pressure ulcer surgery (skin graft or flap). **In this specific case, the use of the low-pressure static mode is preferred.**
- Do not use it with a support with dimensions larger than **120x13x200 cm** or **120x20x200 cm** or **120x25x200 cm**.

# Settings

## 10. Settings

### 10.1 Comfort adjustment

By default, the comfort pressure has been calibrated with a moderate value that is defined on the control panel as **“ STANDARD > 5th LED ”**.

Once the patient is lying on the support, depending on their feelings, their state of comfort and their positioning, you can adapt this pressure using the button  or  PRESSIONS > 5 LED

- i** The pump supervises pressure changes and automatically adjusts them.  
If the pressure drops below the user-defined level, it adjusts the pressure so that it matches the user-set level again. As soon as it is reached, the pump goes into standby mode until the next cycle change.  
If the pressure is too high, it lets the air escape to adjust the pressure.

- i** To obtain an appropriate pressure, refer to the information on page 10 (6.8 manual control) or refer to the indicative table below.

	SOFT		STANDARD		FIRM		
BP <sup>1</sup>	I  II		III  IV		V  VI		SEATED
							
From 30 to 40 kg	From 30 to 65 kg	From 50 to 90 kg	From 60 to 115 kg	From 70 to 140 kg	From 90 to 165 kg <sup>2</sup>	From 120 to 200 kg <sup>2</sup>	From 30° to 90° and 200 kg max <sup>2</sup>

# Settings

The information provided in this table is intended to inform and guide users. They do not constitute a formal recommendation. Users are invited to check the state of comfort as specified in paragraph 6, point 6.8. Nausicaa Medical declines all responsibility for the use or interpretation of the information contained in this table.

**<sup>1</sup> BP = Low Pressure Mode / <sup>2</sup> Depending on the support used!**

## 10.2 Cycle time adjustment

The standard cycle time is calibrated with a moderate value of 15 minutes. The possibility of adjusting the duration can be important depending on the patient's condition, pathology, feeling and comfort, according to medical recommendations or according to the dimensions of the support used.

To change this duration, select the desired value using the (3) button.

Non-exhaustive examples:

- after pressure ulcer surgery > 10 minutes
- support  $\geq 120 \times 13 \times 200$  cm = 20 minutes /  $\geq 120 \times 20 \times 200$  cm = 25 minutes
- Night comfort > 20 or 25 minutes

## 10.3 Mode Setting

### 10.3.1 Dynamic Mode

It is characterized by alternating pressures, thus avoiding prolonged vascular compression that can lead to tissue hypoxia.

As mentioned, (page 33 > 9 Getting Started), this is the mode that is programmed by default. If it is the one that is recommended, you do not have to intervene. It will automatically turn on as soon as the system has reached a minimum pressure of 15mmHg.

**If the pressure of 15 mmHg is not reached within 70 minutes, a visual and audible alarm is triggered to warn of the problem.**

**For more information, see pages 39 & 43.**

### 10.3.2 Static Mode

It allows a greater distribution of the body surface on the support. This is the principle of immersion.

This immersion allows, among other things, the care of people requiring the use of an alternating air support but requiring relative stability of the said support.

This stability is often recommended in the event of immobilization, consolidated fracture, neurological trauma, etc.).

This mode is also required to minimize pain due to local trauma, to promote the patient's rest and during a declared end of life.

# Settings

Depending on the patient's weight, this mode can claim to reduce cell crushing (healing) when the bed bust lift inclination is  $> 30^\circ$ .

For patients who are likely to sleep again on a traditional mattress, whether it is made of standard or memory foam, the static mode can help to carry out withdrawal phases before the installation of the foam mattress.

To turn it on, press the required button (5) to switch from alternating mode to static mode; the control LED lights up to indicate that this mode is on. If necessary, adjust the pressure using the button  or  PRESSION (mmHg)

Press the button (5) again to return to the previous mode.

## 10.3.3 Care Mode

For nursing care, transfer to or from the bed, repositioning, etc., a firm and stable surface of the support is recommended.

To activate this mode, you need to press and hold the button  for more or less 3 seconds. Once activated, it automatically cancels the previous mode, and quickly brings the support to its maximum pressure level, without alternating effect for **30 minutes**.

After this time, to avoid any problems, this function cancels and automatically returns to the level previously defined by the user.

 **For maximum efficiency, it is recommended to activate this mode BEFORE any interventions.**

 It is easy to see if this mode is in operation by looking at the control panel. When this is activated, all LEDs in the comfort range light up one after the other.

 It is possible to disable this mode within 30 minutes. Simply press and hold the button for about 3 seconds.

This action also automatically returns to the previous mode.

 Activating or deactivating this mode is possible even if the control panel is locked.

## 10.3.4 Sitting position

For it to be operational, it is imperative to select this pictogram  before any change in the patient's positioning.

A few seconds after the pump starts, place the patient in a semi-seated position (bed backrest inclined between  $30^\circ$  and  $60^\circ$ ) or sitting ( $60^\circ$  to  $90^\circ$ ). (see photo below).

Here are some recommendations for which this position is very useful:

- To facilitate breathing and blood circulation.
- To allow for better digestion.
- To help prevent pneumonia.
- To keep a horizontal and not vertical visual link.
- To have a meal or swallow a drink while minimizing the risk of choking.

# Settings

• To avoid bottoming out (crushing of cells located in the sacral area) mainly as soon as the bed bust lift reaches an angle  $> 30^\circ$ .

 Dynamic (alternating) mode is typically used with this setting; however, static mode provides a more uniform plan and better comfort.

 **It is important to be aware of the state of safety or danger in certain areas of the body depending on the position a patient occupies.**

- **CERTAIN SAFETY: patient fully lying down ( $0^\circ$ ) up to  $20^\circ$  maximum.**
- **PROBABLE SAFETY: position between  $20^\circ$  and  $30^\circ$ .**
- **DANGER IN SACRED AREA: position between  $30^\circ$  and  $60^\circ$ .**
- **DANGER IN ISCHIAL ZONE: position between  $60^\circ$  and  $90^\circ$ .**

ENGLISH

# Uninstalling

## 11. Uninstalling

- Turn off the pump.
- Remove the plug from the wall outlet and remove the cable from the cable management\* (\*if available).
- Disconnect the support from the pump.
- Remove the pump from the footboard panel.
- Wrap the power cable (e.g. around the housing).
- Replace the pump to its original box (if available).
- Put the **system** back in this carrying bag.

 You can also place the pump in the centre of the support and roll it up. This protects the pump during transport.

# Alarms

## 12. Alarms

This pump has 3 types of alarms.  
These are distinguished by their audible and/or visual signalling.

### 12.1 LOW PRESSURE

It activates to warn that LOW PRESSURE has been detected.  
It is distinguished by an AUDIBLE alarm (a beep at regular intervals) combined with a VISUAL alarm that makes the LED flash at the same rate as the audible alarm.

### 12.2 SERVICE

It manifests itself as soon as the pump requires maintenance.  
It can only be distinguished VISUALLY by means of a fast, consecutive double flashing of the LED at regular intervals.

# Alarms

## 12.3 SYSTEM FAILURE

Indicates severe system failures (software faults or otherwise).

It is characterized by an AUDIBLE alarm (two successive beeps at short and regular intervals) combined with a VISUAL alarm with a double flashing of the LED following the rhythm of the audible alarm.

**⚠ In the event of a maintenance alarm or system failure, the pump must be taken out of service immediately. It is advisable to contact Nausicaa Medical's after-sales service or that of the local service provider for the other steps.**

**i** For hearing comfort, it is possible to deactivate any audible alarm by pressing the reset button.   
Until the problem is resolved, only the LED will continue to flash.  
You can also turn off the visual alarm in the event of a service alarm.

**i** When the problem is resolved, reactivate the alarms by pressing the reset button. 

# Maintenance

## 13. Maintenance

It is **MANDATORY** to clean and disinfect the pump in accordance with applicable national and/or local guidelines for cleaning and disinfecting medical devices.  
In this section, we describe the procedures to be undertaken. It is important to follow them.

**i** To maintain a maximum state of hygiene, the task of cleaning the pump housing is required at the bare minimum of once a week.

**i** When daily or weekly maintenance is carried out on the pump, it is not necessary to disconnect the pump bracket even if it is switched off.

**⚠ The pump is only protected against liquid splashes, as a safety measure during cleaning and disinfection, it must always be disconnected from the mains to avoid any risk of electric shock, fire, burns and consequently, to prevent malfunction.**

**⚠ The pump can carry germs, we recommend that you wear protective gloves BEFORE carrying out any handling, cleaning or disinfection treatment.**

### 13.1 Daily or weekly maintenance

#### 13.1.1 Procedure

- Turn off the pump (mains light off) and leave the support connected.
- Check the pump for external damage.
- Soak a wipe, clean cloth, or single-use paper with a detergent/disinfectant solution and clean the outside of the case. If necessary, let the solution act according to the usual recommendations.

# Maintenance

 **Be careful not to insist too much on the control panel or it will be damaged.**

- Make sure all areas are clean then wipe the case with a clean, dry cloth and allow to dry.
  - Visually check the condition of the filter. It should not have excess dust. If this excess is low (slightly grey filter), remove the filter from its housing and clean with a dust suppressant spray or using low pressure compressed air.
- In the event of more substantial soiling (dark grey or black), it is imperative to provide for its replacement.
- Once the pump has been cleaned, disinfected, dried and inspected, before it is put back into service, make sure that there is no residual moisture on the plug pins, then plug it in and test it to ensure that it is working properly.

## Pre-storage maintenance

### 14. Pre-storage maintenance

#### 14.1 Generality

It is recommended that any defective items be replaced, and that full maintenance be performed before storing the pump. This handling must be carried out by qualified personnel.

#### 14.2 Procedure

 Before any handling and/or intervention, wear protective gloves and carry out an initial surface disinfection.

- Check the integrity of the power cable. It should not show any visible damage or signs of significant wear.
- Check the plug and its pins.
- Check if the case is completely intact.
- If the pump is in order, move it to the cleaning area and place it on a stable, firm work surface.
- Soak a wipe, clean cloth or single-use paper with a detergent solution and clean the entire housing and power cable except the plug.

 **Be careful not to put too much emphasis on the control panel and the electrical pins on the plug, as this may damage them.**

- Make sure the entire case is clean.
- Make the disinfection of the box and electrical cable according to the same process as cleaning.

 **Scrupulously follow the recommendations for the use of the disinfectant solution.**

- Once the afterglow time is over, wipe the remaining excess solution on the two elements (box and electrical cable) with a clean, dry and single use cloth, if necessary, then let dry.
- Regardless of its condition, replace the filter. To do this, remove the filter holder from the back of the pump, remove the used filter, clean and dry the filter holder, insert the new filter and replace the filter holder.

# Pre-storage maintenance

- Once the pump has been cleaned, sanitized and dried, make sure that no residual moisture remains on the pins of the plug.
- Check the internal components for defects and that their fastening system is intact. If necessary, repair and/or replace the defective part(s).
- Connect and test the pump to ensure proper operation.
- After checking that it is working, unplug it, put it back in its box (if available) and place a visible identification label ensuring that it is cleaned, disinfected and that it is working properly.

## DURING THESE INTERVENTIONS, MAKE SURE TO:

- DO NOT submerge or soak the pump.**
- DO NOT spray liquid directly onto the surface of the pump to avoid degradation by liquid penetration inside the housing.**
- DO NOT use a cleaning solution based on hypo carbonate, phenolic, corrosive products, solvents such as acetone, ether, or colouring products (iodine alcohol, etc.) as this could damage the case.**
- DO NOT use abrasive materials.**
- DO NOT leave excess cleaning and/or decontamination solution on the control panel as this may cause damage if penetrated.**
- DO NOT place the pump on the ground.**
- DO NOT wash the filter.**

## Storage

It is recommended:

- Store the pump in a carrying bag or in the original packaging.
- Handle the pump with care.
- Do not put heavy objects on it.
- To protect it from direct sunlight.
- To respect the climatic conditions of storage.

## Climatic requirements

Conditions	Temperatures (°C)	Humidity level (%)	Atmospheric pressure (hPa)
<b>Use</b>	From 10°C to 40°C	From 30% to 75% (non-condensing)	From 700 hPa to 1.060 hPa
<b>Storage / Transport</b>	From 0°C to 40°C		

# Compatible supports

Air height in cm	Dimensions in cm ( L x la )	Foam base	Required cycle time (in minutes)
<b>13</b>	70 x 200	<b>YES</b>	<b>10/15/20/25</b>
	88 x 200		
	120 x 200		
<b>20 (all air)</b>	85 x 200	<b>NO</b>	<b>10/15/20/25</b>
	120 x 200		<b>20/25</b>

# Troubleshooting

## 18. Troubleshooting

The following table lists some problems, possible causes, and the solution to be provided. If the solutions in the panel do not solve the problem, contact your distributor or Nausicaa Medical's after-sales service.

Failure and/or Alarm	Possible cases	Solutions
The pump is not working	The plug is not inserted into the wall outlet or only partially	Push the plug all the way into the wall outlet
	The power outlet is faulty, or the power supply is out of order	Change the outlet or restore the power supply.
	The on/off button is not engaged	Press the on/off button
	The power cable is damaged	Replace the power cable
	1 or 2 out-of-order fuses	Changing the fuse(s)
	<b>If the power light is on and none of the above causes are found, contact your local service provider or Nausicaa Medical's after-sales service.</b>	
No alternating pressure	Static mode is selected	Select Dynamic Mode
	Internal problem (mechanical or other)	<b>Contact your local service provider or Nausicaa Medical's after-sales service</b>
Permanent high pressure	Internal problem (mechanical or other)	<b>Contact your local service provider or Nausicaa Medical's after-sales service</b>
No or little air at the pump outlet	Dirty air filter Compressor out of order.	Check the condition of the air filter Changing the compressor
Service Alarm	The pump requires maintenance	<b>Contact your local service provider or Nausicaa Medical's after-sales service</b>
System failure alarm	Internal problem (mechanical or other)	<b>Contact your local service provider or Nausicaa Medical's after-sales service</b>
Low pressure in the cells of the support hence Low-pressure alarm	- The set pressure is too low in relation to the patient's weight and/or position (seed-sitting or sitting). - The patient's weight is greater than that allowed by the pump	- Increase the pressure - Adapt according to position - Use a pump that supports the patient's weight
Low Pressure Alarm	CPR open Air leakage at connections Cells leakage	Close the CPR Check Connections Check Cells & Connections

# Symbols (product label and/or packaging)

## 19. Symbols (product label and/or packaging)

	Indicates that the instructions for use should be read		Indicates an important warning
	Advice or information that the user and/or patient should be aware of		Indicates that the instructions for use should be read before using the product
	Indicates that the product is a medical device that complies with the EU Medical Device Regulation 2017/745		Indicates that the product complies with European Union legislation
	Unique identification code for medical devices		Manufacturer's Contact Information
	Date of manufacture		Serial Number
	Article reference		Responds to the norm CEI 60601-1 - Type BF
<b>IP21</b>	Ingress protection		Indicates fuse type
	Electrical insulation class II: double insulation		Do not dispose of household waste
	Do not use chlorinated agents (Cl)		Patient weight

# Technical data

## 20. Technical data

Model	R'GO SOINS
Article reference	RGOSOINS-COMP-PM
MDR Class	Class I
IUD basic	37014294SYSRGOSOINSTX
UDI-DI	03701429409050
Electrical protection class	 Class II, Type BF, (No AP/APG)
Noise pollution	27-28 dB (Maximum flow rate at 1 m) (ISO 7779)
Mains voltage	230 V
Frequency of the sector	50 Hz
Fuse Rating	250 V, T1A
Ingress protection	IP21
Airflow	7 L/min.
Pressure Range	15 to 50 mmHg
Cycle times	10/15/20/25 minutes
Housing material	ABS plastic
Dimensions	(L) 280 mm x (L) 205 mm x (H) 135 mm
Weight	2,6 kg

## Elimination

### 21. Elimination

Comply with national regulations to remove any of the components of this pump or to eliminate it completely.

## Lifespan

### 22. Lifespan

The lifespan shown below depends on the frequency and intensity of use, cleaning, disinfection, transport, compliance with annual maintenance, etc.

Product	Lifespan
R'GO SOINS pump	± five (5) years

**The fact that Nausicaa Medical specifies a lifespan does not constitute an additional guarantee.**

# Warranty

## 23. Warranty

**Article 1 :** NAUSICAA Medical S.A. warrants this device against any manufacturing defects and assembly of its components, and only for devices used under the conditions provided by NAUSICAA Medical S.A. This warranty, the conditions of which are defined below, is valid for 24 months from the date of the first purchase from NAUSICAA Medical S.A.

**Article 2 :** The warranty entitles you to free labour and free replacement of parts found to be defective.

**Article 3 :** The outbound port of the device, as well as all related costs, are the responsibility of the reseller. The goods always travel at the risk and under the responsibility of the dealer.

Sous garantie : the costs of return after intervention will be borne by the company NAUSICAA Medical S.A.

Hors garantie : the return costs are at the expense of the retailer whether or not he accepts the repair estimate.

**Article 4 :** The warranty does not apply if the claims are resulting from:

- Accident, misuse of the appliance, or negligence of the purchaser.
- Carrying the device without adequate protection.
- A modification or transformation not validated by the company NAUSICAA Medical S.A.
- The impact of external agents (natural disaster, fire, shocks, humidity, floods, lightning, etc.).
- Installation and/or use in a manner that does not comply with technical and safety standards if the device is to operate in a country other than the country of purchase and/or if the power supply is not suitable for the operating voltage of the device.
- A lack of routine maintenance.

**Article 5 :** The reseller will not be able to invoke the benefit of the warranty:

- If the device's serial number has been removed, altered, or rendered illegible.
- If the device under warranty has been modified without the approval of NAUSICAA Medical S.A.

**Article 6 :** During the repair of defective equipment, no equipment loan will be made.

**Article 7 :** Any warranty claim shall be exercised through the reseller. If this is not possible, the buyer may send his equipment directly to the company NAUSICAA Medical S.A. In this case, the buyer must indicate in a letter attached to the equipment, the contact details of the dealer and a copy of the purchase invoice.

**Article 8 :** The shipment of spare parts under warranty will be made or will only be made after consultation with the After-Sales Service of NAUSICAA Medical S.A.

Please note that defective spare parts must be returned to the After-Sales Service of NAUSICAA Medical S.A. within 15 days, otherwise they will be invoiced 1 month after the parts have been sent by NAUSICAA Medical S.A.

**Article 9 :** Defective parts changed under or out of warranty are guaranteed for 6 months from the date of repair or shipment of spare parts.

**Article 10 :** No dealer may unilaterally change the terms of this warranty.

# NAUSICAA

MEDICAL

Sales Department Headquarters

Phone : 04 66 51 50 80

Fax : 04 66 51 50 47

Mail : [contact@nausicaa-medical.com](mailto:contact@nausicaa-medical.com)

[www.nausicaa-medical.com](http://www.nausicaa-medical.com)

Your contact :