

PATIENT LIFT FLYER 300

With electrical opening legs

Nausicaa Médical offers you a range of XXL products. With a maximum capacity of 300 kg, all our bariatric products can be used for transfers and to care for bedridden patients. The FLYER 300 is a comfortable, yet powerful lift, with an XXL U-shaped sling for rapid positioning by the caregiver.

Easy to use, it provides safe, comfortable transfers for patients.

It comes with a lithium battery as standard, avoiding additional operating costs.

Ses + sont :

- Easy to drive
- Comfort flail
- Very wide range of movement
- Electric foot spreader
- Reinforced chassis for bariatric patients


Features



CONFORT STYLE BIEN-ÊTRE
FABRICATION FRANÇAISE

Dimensions

- A. Total length: 134 cm
- B. Minimum useful height: 53,5 cm
- C. Max. usable height: 170,5 cm
- D. Total height: 127 cm
- E. Frame height: 11,5 cm
- F. Minimum base width: 78 cm
- G. Max. baseplate width: 101 cm
- H. Turning diameter: 142 cm



REFERENCE

FLYER300-CBLI