

U-SHAPED SLING XXL

The positioning of the shoulder hooks, which run along the entire back of the sling, prevents compression of the patient during transfer.

Standard comfort foam in the leg slings reduces shearing.

The U-shaped design prevents the slings from being placed under the patient's buttocks, making donning and doffing more comfortable for the patient and easier for the caregiver. This sling must not remain under the patient.

- Available only in comfort jersey
- Sizes 2XL, 3XL, 4XL.

ADJUSTING HANDERS



The multi-hanger system allows you to vary the patient's position. For optimum sitting comfort, we recommend setting the sling as short as possible on the shoulders (yellow setting) and as long as possible on the legs (black setting).

SLINGS POSITIONING



CROSSED BETWEEN LEGS

Suggested choice



CROSSED UNDER THE LEGS

Secondary choice

HEADREST (OPTIONAL)

Optional headrest with fasteners (TAAT-PF) and removable headrest (TAMOV-PF).



CONFORT STYLE BIEN-ÊTRE
FABRICATION FRANÇAISE



REF.	SUCBAS-PF
Option	TAAT-PF (headrest with fasteners) TAMOV-PF (removable headrest)